

14 STRONGARM™ SURGICAL HOLDERS AND POSITIONERS

TWIN ARM SERIES
ONE DEVICE - TWO FUNCTIONS



■ **REF. NO. 73010**
Twin Arm StrongArm Holder and Positioner with Quick Disconnect feature on both arms and complete with a Heavy Duty Rail Clamp.

■ **REF. NO. 73012**
Twin Arm StrongArm Holder and Positioner with Quick Disconnect feature on both arms. (No Rail Clamp.)



■ **REF. NO. 73040**
StrongArm-Flex Arm Combo System with Quick Disconnect features on both arms and complete with a Heavy Duty Rail Clamp.

Dimensions for Twin Arm Units (73010 and 73012):

Post Height: (A to B) Adjustable from 3.5" - 11.5" (9cm - 29cm)
Mid Segment (Both Arms, B to C): 10.5" (26.5cm)
End Segments (Both Arms, C to D): 13" (33cm)

Note: Tips added to the Quick Disconnect Fitting Add about 3" (7.5cm) to the End Segment Length

Dimensions for Combo System:

Same as above except for FlexArm dimensions (See page 5):

THE STRONGARM-FLEX ARM COMBO SYSTEM IN USE.

The StrongArm Portion is holding an available Mediflex Nathanson Hook (See pg 16) while the Flex Arm Portion is Holding a 10 mm Laparoscope (page 27) with the available 10mm Quick Grip Scope Holder (page 15). This features the near absolute rigidity of a StrongArm Device with the flexible rigidity of a Flex Arm Holder and Positioner. This combination is ideal for procedures that require the liver to be retracted out of the way while also holding a scope. It can allow the surgery to be done with fewer employees... a true labor saving device.



See pages 15 and 16 for available tips and retractors to fit all Mediflex Holders and Positioners equipped with the Quick Disconnect Feature. See page 15 for available Sterilizing Trays.

